



SurePower™ Battery Pack Compatibility Insert

The ZOLL SurePower™ Battery Pack is compatible with the following ZOLL Defibrillators:

R Series™ Defibrillators

Compatible with *all* R Series Defibrillators.

E Series® Defibrillators

Compatible with E Series® Defibrillators manufactured after **July 2006**, with serial numbers in the format of **AB06Hnnnnnn**, where nnnnnn is a number specific to the defibrillator. E Series Defibrillators must run **Software Version 3.00.00** or higher to be compatible.

E Series Defibrillators manufactured before July 2006 may be returned to ZOLL Medical Corporation to be upgraded for use with the SurePower Battery Pack.

AED Pro® Defibrillators

Compatible with AED Pro® Defibrillators with serial numbers higher than AA06E002653. AED Pro Defibrillators must run software version 2.40 or higher to use the SurePower Battery Pack.

AED Pro Defibrillators with serial numbers lower than AA06E002653 can be returned to ZOLL Medical Corporation and upgraded for use with the SurePower Battery Pack.